



**Buy 1 case of KIMBLE® GLS 80® Bottles  
Receive a FREE pack of Four Port GLS 80 caps**

### KIMBLE® GLS 80® Wide Mouth Media Bottles w/ Blue PP Cap

KIMBLE® GLS 80® wide mouth media bottles are extremely versatile: ideal for pouring, filling, mixing, preparing, storing, transporting, dispatching, taking samples in the field and hot filling, whether in the laboratory or in production. The wide mouth design facilitates easy and safe charging and removal of powders, pastes, granules and liquids. Spoons, spatulas and tweezers can be manipulated with ease. Even large-volume funnels can be inserted into the bottles without any difficulty. The wide neck also makes cleaning the bottles particularly easy.

Durable, fired on ceramic ink graduations and marking patch withstand the rigors of daily laboratory, bioprocess and production use. KIMAX® 33 expansion borosilicate glass is low extractable, USP Type I, your assurance of purity, with a full temperature and chemical resistance range. Complete with blue quick-release linerless polypropylene closure and pour ring that are rated up to 140°C. The exterior tint on amber bottles provides UV light protection in the 300 nm – 500 nm wavelength range. Autoclavable KimCote® coating provides safety and reduces hazards of scratching, leaking and shattering.

Clear Cat. No.	Amber Cat. No.	KimCote® Cat. No.	Capacity (mL)	OD (mm)	H (mm)	Qty/Case
14393-0250	14394-0250	14390-0250	250	95	105	4
14393-0500	14394-0500	14390-0500	500	101	152	4
14393-1000	14394-1000	14390-1000	1000	101	222	4
14393-2000	14394-2000	14390-2000	2000	136	252	4
14393-5000	14394-5000	14390-5000	5000	182	314	1
14393-10L	14394-10L	—	10000	227	389	1
14393-20L	14394-20L	—	20000	288	484	1



### GLS 80® Four Port Screw Cap\*

The multi-port system consists of a GLS 80 screw cap with four ports. The tubes\* are fed through these ports and secured with inserts\* and screw connections\*. There is a choice of 5 inserts, each with different internal diameter depending on the application. Unused ports can be fitted with a blanking cap, if necessary. Sterile pressure compensation can be achieved using the membrane filter set (0.2 µm pore size).

Cat. No.	Thread GLS	d mm	h mm	Qty/Case
293109103	80	86	~60	2
<b>*Accessories (sold separately)</b>				
1160169	Insert for screw cap GL 18, ID 3.2mm			1
1160170	Insert for screw cap GL 18, ID 6.0mm			1
1160171	Insert for screw cap GL 18, ID 8.0mm			1
1160172	Insert for screw cap GL 18, ID 10.0mm			1
1160173	Insert for screw cap GL 18, ID 12.0mm			1
1160174	Open-top screw cap for tube connection, blue, GL 18			2
1160167	Pressure equalization kit for 4-port, GL 18			1
1170682	Closed-top screw cap, red, PBT, GL 18			2

\*Tubing, inserts and screw connections must be purchased separately. Visit [www.kimble-chase.com](http://www.kimble-chase.com) to view accessories.





Excellence in your hands

### To claim your redemption (redeem by 1/31/22)

1. Obtain your Proof of Purchase (dated between 7/1/21- 12/31/21). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit [www.dwkdeals.com](http://www.dwkdeals.com)
3. Attach your Proof of Purchase and fill-out the online form.
4. The value of freeware will be calculated off of DWK list price and units of measure found on [www.DWK.com](http://www.DWK.com). Free items must be of equal or lesser value. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. DWK distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.

Contact us at [DWK.com](http://DWK.com)